



Safety Syringe, three-component LuerLock:

The safety syringe from MedicTinedic ApS provides patients as well as personal with an increased level of security and insures correct medical procedures.

The patented pullout system ensures the needle is fixated within the needle after usage. This drastically reduces the risk of accidental needle stick injuries with the subsequent anxiety, examinations and possible leave of absence.

The safety syringe can be delivered as an end-user product with the needle already attached or as a syringe without needle for customized options. The patented system is perfectly functional with all types of needles.

Sizes:

2-3ml

5-6ml

10-12ml

Syringes:

18-27 G

Advantages:

- A disposable syringe that truly ensures single usage.
- After injection the piston is retracted and the needle is fixated with the syringe. For added safety measures the piston can be broken in two and inserted back into the syringe.
- Reduced risk of HIV or hepatitis infection by accidental needle stick injuries.
- Prevents the reuse of syringes and drastically reduces spread of severe illnesses.
- Reduces leave of absent for personal and staff.
- Provides more comfortable working conditions for personal and staff.
- Comparable syringes are equal in price or with neglectable lower costs.

